



medical **ECONET**
GERMANY



ECOspin III (IVDR)

Compact benchtop centrifuge for daily routine

- 12-fold angle motor up to 15ml tubes
- Brushless motor with electronic speed control in steps of 100 rpm
- Digital preselection of running time and speed
- Hygienic membrane keypad
- Electromagnetic lid lock
- Durable, easy-to-clean plastic housing
- In vitro diagnostics (IVDR) certified centrifuge

Specifications

Equipped with a 12-position angle rotor for holding tubes with a volume of up to 15 ml and a max. RCF of 2,400 g, the ECOspin III is the ideal device for general practitioners.

The ECOspin III centrifuge has a small footprint and is easy to operate. The last selected settings for speed and time are memorised. This makes it possible to start at the touch of a button. For the user's safety, the lid is locked electromagnetically during operation. The hard-wearing plastic housing and hygienic membrane keypad make cleaning easy.

Translated with www.DeepL.com/Translator (free version)



Dimensions

255 mm (W) x 245 mm (D) x 140 mm (H)



Specifications

| | |
|------------------------------|-------------------------------------|
| Rotor | 12-place rotor for up to 15ml tubes |
| Speed Control | elect. in steps of 100 rpm |
| Max. RCF | 2.400 g |
| Max. Rotation Speed | 4500 U / min |
| Min. Rotation Speed | 300 U / min |
| Power | 110-240V 50/60 Hz |
| Noise at max. Rotation Speed | approx. 56 dB (A) |
| Dimensions (W x D x H) | 255 x 245 x 140mm |
| Weight | 6 kg |

medical **ECONET**

Im Erlengrund 16, 46149 Oberhausen - Germany

t: +49 (0) 208 377 890-0 info@medical-econet.com
f: +49 (0) 208 377 890-55 www.medical-econet.com
facebook.com/medicaleconetgermany
twitter.com/medicaleconet

Distributor

*Subject to change